# A transformational one-step library preparation method for multiplexed plasmid and amplicon sequencing



Jack T. Leonard, Rebecca Feeley, Stella Huang, Li Wang and Jessica Smith seqWell, Inc., Beverly, MA USA

#### Introduction

Conventional library preparation methods demand multiple incubations interrupted by cumbersome pipetting steps. A transformational new library prep technology from seqWell enables barcoding and amplification of 8 - 384 dual-indexed libraries in 90 minutes in 16 µl "one-pot" reactions, while simultaneously auto-normalizing library read-count and insert size. The simplicity of the one-step workflow makes it uniquely well-suited for manual, automated <u>and</u> ultra high throughput (UHT) library prep from plasmids and PCR products. To demonstrate miniaturization, we prepared plasmid libraries on the SPT LabTech mosquito<sup>®</sup> HV at 3,200 nanoliters per reaction in 384-well PCR plates, and the sequencing data were used for *de novo*, full plasmid assembly.

## One-step library prep workflows





#### B. UHT one-step 384-plex plasmid libraries with variable DNA input: Typical depth of coverage profiles and de novo plasmid assembly success rate



Summary: 382 out of 384 plasmids resulted in circularized de novo assemblies that were 100% concordant with their known references (i.e., 99.5% overall success rate).

### Conclusions

- Full-volume 96-plex one-step reactions (16 µL) set-up manually or with standard automated liquid handlers produce a uniform number of sequencing reads per plasmid (c.v. 16.6%).
- 384 miniaturized one-step reactions (3,200 nL) set-up with an 8-fold DNA input range (1.6 12.8 ng per reaction) resulted in a very high *de novo* plasmid assembly success rate (99.5%), even when sequenced to modest depth (3.3 5.6 K read pairs per sample) on a 2 x 151 cycle MiSeq Micro v2 kit.
- These results suggest that the *de novo* plasmid assembly success rate would be in the vicinity of 99.5%, if 768 one-step plasmid libraries (8 x 96-well plates) were pooled and sequenced on a 300 cycle MiSeq Micro v2 kit, or, if 2,304 one-step plasmid libraries (6 x 384-well plates) were pooled and sequenced on a 300 cycle MiSeq V2 kit.