



MosaiX™ Library Prep Kit

Simplicity of tagmentation,
complexity of ligation


High-Performance TnX Library Prep




MosaiX Applications



Whole-Genome Sequencing



Low-Pass Whole-Genome Sequencing



Metagenomic Sequencing



Whole-Exome Sequencing



Targeted Hybrid Capture Sequencing

No need to choose between ease-of-use, library complexity or cost-efficiency. With MosaiX you can have all three.



High-Performance

- TnX next-generation transposase
- Uniform sequence coverage
- Increased library complexity
- Low duplication rates



Rapid

- 90-minute workflow
- Eliminates DNA shearing
- Automation-friendly



Flexible


- Indexing strategy
- Batch size & DNA input
- Buffer tolerance



Cost Savings

- Affordable per sample pricing
- Reduced time & labor
- Fewer consumables

seqWell Directional Tagmentation with TnXcomplexity made simple



Better transposase Better libraries Better NGS

Generate high-quality libraries with TnX next-generation transposase – engineered for reduced insertion bias, increased transposase activity, and enhanced inhibitor tolerance to provide robust performance and improved uniformity of coverage.




Figure 1: Read start base bias for bead-linked tagmentation (Tn5), Mosaix (TnX), and enzymatic fragmentation are shown. Reduced bias of TnX versus Tn5 is observed with a profile similar to that of enzymatic fragmentation.




Figure 2: Mosaix fragmentation profiles are unaffected by EDTA allowing DNA dilution in a variety of buffers.

seqWell Directional Tagmentation with TnXcomplexity made simple

MosaiX library prep generates high-complexity libraries with an automation-friendly workflow powered by TnX, the next-generation transposase engineered for reduced insertion bias, increased enzyme activity, and enhanced inhibitor tolerance. TnX simplifies fragmentation by eliminating the need for time-consuming DNA shearing or digestion while providing improved uniformity of coverage compared to alternative library prep methods.

Streamlined, 90-minute Workflow




Figure 3: Times for technology workflows are based on manufacturers' standard protocols. During seqWell directional tagmentation, the tagging and ligation reactions occur sequentially within the same well. Following bead clean up, P5 and P7 indexes are added via PCR. **Inset** – Following PCR, libraries prepared from 1-50 ng of human genomic DNA (NA12878) using MosaiX were diluted 1:10 and run on a TapeStation D5000 HS kit for fragment sizing. Mean fragment sizes (with adapter) of ~550 – 800 bp were seen independent of DNA starting concentration.




"I have been really impressed with the **outstanding performance** of MosaiX library prep... ..(which) **reliably resulted** in libraries with **high complexity and uniform coverage** that provided for more **precise sequencing data**. All this was achieved using an extremely **efficient workflow**."

Massimo Delledonne, PhD

Professor of Genetics, University of Verona, Italy

High-quality Libraries Using 50% Shorter Workflow

Superior Performance & Scalability - Whole Genome Sequencing (WGS)




MosaiX library prep led to superior WGS sequence data compared to a competitor tagmentation method, including higher mean coverage, lower duplication rate and larger estimated library size (complexity). Its streamlined workflow provides the speed and scalability to support WGS efforts of any size.

Library Prep Method	Mean Coverage (X)	% Bases $\geq 20X$	Duplication (%)	Estimated Library Size
MosaiX	27.4	84%	10%	3,374,667,137
Bead-linked Tn5	23.9	80%	14%	2,871,092,622

Figure 4: Libraries were generated from 40 ng of NA12878 DNA (Genome in a Bottle) using MosaiX (seqWell) or a competitor tagmentation-based kit following manufacturers’ standard protocols. Libraries were sequenced on a NovaSeq™ X Plus, data down sampled to 105 Gb, and sequences aligned to hg38. A normalized coverage plot including mean base quality for MosaiX libraries shows the reduced GC bias of TnX.

Find Those Critical Exome Targets with TnX

Improved Variant Detection - Whole Exome Sequencing (WES)



No need to sacrifice variant calling quality to scale up your WES projects. TnX transposase can access difficult regions that Tn5 cannot, leading to less target drop outs compared to competitor transposase methods.

Library Prep Method	Duplication (%)	HS Library Size	% Zero Cvg
MosaiX	4.2%	285,255,870	0.69%
Bead-linked Tn5	7.9%	167,531,020	0.82%
Enz. Frag & Ligation	6.0%	201,223,608	0.68%

Figure 5: Libraries were generated using MosaiX, competitor bead-linked tagmentation, or competitor enzymatic fragmentation and ligation from 40 ng of NA12878 DNA (Genome in a Bottle). Libraries underwent target capture using Twist’s Exome 2.0 panel and sequenced on a NextSeq™ 2000. Data were down sampled to 6 Gb each and aligned to Twist exome capture targets on hg38.

Specifications

Sample Input Types	Genomic DNA
DNA Input Range	1 - 40 ng*
Buffer Compatibility	DNA diluted in 10mM Tris-HCl, 1X TE, Low TE, or Water
Total Library Prep Time	90 minutes (35 minutes hands-on time)
Mean Output Fragment Size** <i>Fragment sizes as measured by TapeStation (including adapters)</i>	<ul style="list-style-type: none"> • WGS protocol (recommended 0.75x bead clean up): 700 bp \pm 15% • Target capture protocol (recommended 0.8x bead clean up): 650 bp \pm 15%
Indexing Strategy	Any tagmentation-compatible indexing or custom primers can be used
Supported Paired Reads (clusters/sample)	\geq 400 million
Reactions per kit	24 or 96
Sequencer Compatibility	<ul style="list-style-type: none"> • All Illumina sequencing platforms • Compatible with Complete Genomics and BGI platforms or Element Biosciences AVITI™ using conversion kits for Illumina libraries

* <5 ng may require optimization of adapter concentration and PCR cycles

Ordering Information

Catalog Number	Description	Number of Reactions
301458	MosaiX Library Preparation Kit	24 reactions
301464	MosaiX Library Preparation Kit	96 reactions
301443	MosaiX Compatible UDI Primers - Set 1	96 primers
Inquire	For purchase of more than 96 primers	

www.seqwell.com/mosaix
sales@seqwell.com

Request a Quote



© 2026 seqWell, Inc. All rights reserved.
 For Research Use Only. seqWell, the seqWell logo, TnX, and MosaiX are trademarks of seqWell, Inc. All other trademarks are property of their respective owners.